

# PHILIPS

## MR Patient Care

# Philips MR Patient Monitoring solutions

### Premium MR Patient Monitor for high-acuity care Expression MR400



- 15" LED touchscreen
- 5,000 Gauss
- 4 W/kg SAR & 7.2  $\mu\text{T}$  B1<sub>rms</sub>
- 1.5T & 3.0T MR compliance
- Gauss alarm (with AA option)
- Parameters:
  - Basic: ECG\*, SpO<sub>2</sub>\*, NIBP, CO<sub>2</sub>/Respiration
  - Optional: IBP (x2), body temperature, anesthetic agents
- Advanced clinical alarms
- MR gating feature: wired/wireless

\* wireless ECG & SpO<sub>2</sub>

### Premium MR Patient Monitor for low-acuity care Expression MR200



- 12" LCD
- 1,500 Gauss
- 4 W/kg SAR
- 1.5T & 3.0T MR compliance
- Gauss alarm (standard)
- Parameters:
  - Basic: ECG\*, SpO<sub>2</sub>\*, NIBP
  - Optional: CO<sub>2</sub>/respiration
- MR gating feature: wired/wireless

\* wireless ECG & SpO<sub>2</sub>

### Information Portal

### MR Patient Care Portal 5000



- Medical grade projected capacity touch panel 18.3" LED display
- Remote monitoring with Expression MR200 and Expression MR400
- Modular design for cost effective and possible future upgrades
- HIS/EMR connectivity with enhanced product and IT security
- Patient management
- Password protected administrative profile features
- Advanced visual perception principles and decision-making tools
- Optional: barcode scanner, printer, wireless keyboard and mouse

**More than 30 years of excellence in MR patient monitoring**

### Key facts at a glance

#### Patient types

- Adults
- Infants
- Pediatric
- Neonates



#### Parameters & alarms

- Bedside-level
- Advanced clinical alarms



#### Make the most of the MR

4 W/kg SAR | 5,000 Gauss\*  
7.2  $\mu\text{T}$  B1<sub>rms</sub>\* | 1.5T & 3.0T

#### Connect

- with convenience
- Enhanced product and IT security
- Wireless communication
- HL7/RS232 connectivity HIS/EMR



\* Based on Expression MR400 specifications

Your Philips contact:

